



Be Right™



Conductivity Standard Solution, 1000 $\mu\text{S}/\text{cm}$ (NIST), 491 mg/L NaCl, 100 mL

Product #:

1440042

AED Price:
Available

Contact Hach

491 mg/L NaCl, for conductivity probes,

also suitable for standardising conductivity apparatus, 1000 $\mu\text{S}/\text{cm}$ ($\mu\text{mhos}/\text{cm}$)

Specifications

Concentration:	491 mg/L NaCl
Parameter:	Conductivity
Platform:	Solution
Quantity:	100 mL