



NitraVer 5 Nitrate Reagent Powder Pillows, 0.1-10mg/L NO₃-N

Product #:

1403499

AED Price:
Available

Contact Hach

Hazardous



Powder reagent in sealed packages for the determination of nitrate. For 25 mL sample, 100 pcs.

Specifications

Description:	NitraVer 5 Nitrate Reagent Powder Pillows, 25mL
Instrument:	DR/3000, DR/2000
Measuring range:	0.3 - 30.0 NO ₃ -N
Measuring range (2):	0.1 - 10.0 mg/L NO ₃ -N
Method:	Cadmium Reduction
Method number:	8039 HR
Method number :	8171 MR
Number of tests:	100
Parameter:	Nitrate
Platform:	Powder Pillows
Quantity:	100 /pk
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)