



NitriVer 3 Nitrite Reagent Powder Pillows

Product #:

2107169

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Nitrogen, Nitrite 0.002-0.300mg/L NO₂⁻-N. For 10 mL sample, 100 pieces.

Specifications

Description:	NitriVer 3 Nitrite Reagent, 10 mL
EPA compliant:	Yes
Instrument:	All except PC II
Measuring range:	0.002 - 0.300 mg/L NO ₂ ⁻ -N
Method:	Diazotisation
Method number:	8507
Number of tests:	100
Parameter:	Nitrite
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)