



Nitrogen, Ammonia TNT, 0.02 - 2.50 mg/L NH₃-N

Product #:

2604545

Hazardous

AED Price:
Available

Contact Hach



Cuvette test (Ø 16 mm) completely equipped with all reagents, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube. 25-50 tests

Specifications

Description:	Nitrogen, Ammonia Reagent Set, AmVer (Salicylate), LR
Instrument:	DR 900
Measuring range:	0.02 - 2.50 mg/L NH ₃ -N
Method:	Salicylate
Method number:	10023
Number of tests:	25 - 50
Parameter:	Ammonia
Platform:	Test 'N Tube
Special remarks:	Set contains: Cuvette test tubes, AmVer Diluent Reagent (2602200), Ammoniasalicylate Powder Pillows (2395266), Ammoniacyanurate Powder Pillows (2395466)
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)