



Nitrite TNT, 0.002 - 0.500 mg/L NO₂-N

Product #:

2608345

Hazardous

AED Price:
Available

Contact Hach



Cuvette test (Ø 16 mm) completely equipped with all reagents, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube
For determination of nitrite by diazotization. 50 tests

Specifications

Description:	NitriVer 3 TNT Reagent Set
Instrument:	DR 900
Measuring range:	0.003 - 0.500 mg/L NO ₂ -N
Method:	Diazotization
Method number:	10019
Number of tests:	50
Parameter:	Nitrite
Platform:	Test 'N Tube

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)