



Nitrogen, total, TNT, 0.5-25 mg/L N

Product #:

2672245

Hazardous

AED Price:

Contact Hach

Ships within 3 days



Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Description:	Nitrogen, total, TNT Reagent Set
Digestion required:	Yes
Instrument:	DR 900
Measuring range:	0.5 - 25.0 mg/L N
Method:	Persulphate Digestion
Method number:	10071
Number of tests:	25 - 50
Parameter:	Nitrogen, total
Platform:	Test 'N Tube
Special remarks:	Set consists of 2 cartons: 26717-45 and 26721-45 26717-45 = TN Hydroxide reagent (2714045), Persulfate Powder Pillows (2671846) 26721-45 = TN Acid reagent C (2672145), Bisulfite reagent A Powder Pillows (2671946), Indicator reagent B Powder Pillows (2672046)
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)