



## Ozone Accuvac, LR, 0.01-0.25 mg/L O<sub>3</sub>, pk/25

**Product #:**

**2516025**

**Hazardous**

AED Price:  
Available

Contact Hach



Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of ozone by the Indigo method. pk/25

---

### Specifications

Description:	Ozone AccuVac, LR
Instrument:	All except DR/820
Measuring range:	0.01 - 0.25 mg/L O <sub>3</sub>
Method:	Indigo
Method number:	8311
Number of tests:	25
Parameter:	Ozone
Platform:	Accuvac

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)