



Phosphate, total, TNT, 1.0 - 100 mg/L PO₄

Product #:

2767245

Hazardous

AED Price:
Available

Contact Hach



Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Description:	Phosphate, total,TNT Reagent set, HR
Digestion required:	Yes
Instrument:	DR 900
Measuring range:	1.0 - 100 mg/L PO ₄
Method:	Molybdovanadate with Acid Persulfate Digestion
Method number:	10127
Number of tests:	25 - 50
Parameter:	Total Phosphorous
Platform:	Test 'N Tube
Special remarks:	Set contains: Cuvette test tubes, Potassiumpersulfate Powder Pillows (2084766), Molybdovanadate Reagent (2076026), Sodiumhydroxide solution 1.54 N (2743042), Sulfuric acid solution 1.37 N (2742900)
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)