



Molybdovanadate Reagent Solution (100 mL MDB)

Product #:

2076032

Hazardous

AED Price:
Available

Contact Hach



Molybdovanadate Reagent for the determination of ortho Phosphate, 0.3 - 45.0 mg/L PO₄

Specifications

Instrument:	All except DR/850, DR/820, PC II
Measuring range:	0.3 - 45.0 mg/L PO ₄
Method:	Molybdovanadate
Method number:	8114
Number of tests:	50
Parameter:	ortho Phosphate
Platform:	Solution
Quantity:	100 mL
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)