



Phosphate Reagent Powder Pillows 0.02-2.50mg/L PO₄

Product #:

2106069

Hazardous

AED Price:
Available

Contact Hach



Powder reagents in sealed packages for the determination of Ortho-Phosphate. For 10mL samples; pk/100.

Specifications

EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000, DR/820
Measuring range:	0.02 - 2.50 mg/L PO ₄
Method:	Ascorbic Acid
Method number:	8048
Number of tests:	100
Parameter:	Phosphate, ortho
Platform:	Powder Pillows
Quantity:	100 /pk
Reference method:	SM 4500-P-E
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)