



PhosVer 3 Phosphate Reagent Powder Pillows 0.02-2.50mg/L PO₄

Product #:

212599

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Phosphorus, reactive (25mL sample; 100 tests).

Specifications

Instrument:	Best for use with DR2000, DR3000, DR900; or on other spectrophotometer with use of the 25-mL sample cell.
Measuring range:	0.02 - 2.50 PO ₄
Method:	Ascorbic Acid
Method number:	8048
Number of tests:	100
Parameter:	Phosphate
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)