



Maintenance kit for 5500 sc Silica Analyser, 2/4 channel

Product #:	6788302
AED Price:	Contact Hach

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Specifications

Product Kit:	5500 sc Silica Analyser, 2/4 Channels
What's included?:	Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

What's included?

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions