



## Maintenance kit for 5500sc Silica Analyser, 6 channels

Product #:	6788303
AED Price:	Contact Hach

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

---

### Specifications

Product Kit:	5500 sc Silica Analyser, 6 Channels
What's included?:	Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

---

### What's included?

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions