



## Pipette, variable, volume 1.0-5.0 mL

<b>Product #:</b>	<b>BBP065</b>
AED Price:	Contact Hach
Available	

Hach's variable volume single-channel, air-displacement pipettes feature robust design, comfortable ergonomics, slim-line tip ejection, and single-handed operation. Outstanding accuracy with adjustments in 50 µL increments (0.05 mL). Ideal for use with Hach LCK reagents. Supplied with calibration report and maintenance tool.

---

### Specifications

Parameter:	Pipette
Quantity:	1 each
What's included?:	Pipette 1.0-5.0 mL, maintenance tool, calibration certificate, basic manual.  Tips are not included and have to be purchased separately.

---

### What's included?

Pipette 1.0-5.0 mL, maintenance tool, calibration certificate, basic manual. Tips are not included and have to be purchased separately.