



## Pocket Pro+ pH Tester with Replaceable Sensor

**Product #:** 9532000  
**AED Price:** Contact Hach

### Take The Guesswork Out of Your Measurements.

Hach Pocket Pro and Pro+ are engineered to deliver accurate results. Backed up with built in performance diagnostics, you never have to guess when to clean or calibrate the sensor. You also get replaceable batteries for convenient field use, and a large, easy-to-read LCD screen.

#### Confidence in results

Easy calibration steps and built-in diagnostics for pH testers take the guesswork out of calibrating the sensor to keep it in optimum condition.

#### Saves time and money

Easy-to-replace AAA batteries and a convenient sample cup make maintenance and measurements simple.

#### Ease of use

The large LCD screen is easy to read; backlighting is included on Pocket Pro+ models, allowing measurement in almost any lighting environment.

---

### Specifications

Accuracy:	± 0.02 pH
Backlight:	Yes
Battery requirements:	4 AAA
Casing IP rating:	IP67
Manual languages:	English, French, German, Italian, Spanish
Operating temperature range:	0 - 50 °C
Parameter:	pH, Temp.
Power Supply:	Batteries (see Battery Requirements)
Range:	0.00 - 14.00 pH
Resolution:	0.01 pH
Sample volume:	9 mL
Warranty:	24 months
What's included?:	Includes: tester, 4 AAA batteries, lanyard and quick start guide.

---

### What's included?

Includes: tester, 4 AAA batteries, lanyard and quick start guide.