



Pocket Pro+ ORP Tester with Replaceable Sensor

Product #: 9532100
AED Price: Contact Hach

Take The Guesswork Out of Your Measurements.

Hach Pocket Pro and Pro+ deliver accurate results with built in diagnostics, replaceable batteries, and a large, easy-to-read LCD screen.

Confidence in results

Take the guesswork out of when to calibrate or clean the sensor with built-in diagnostics

Saves time and money

Easy to replace AAA batteries and convenient sample cup make maintenance and measurements simple

Ease of use

Large LCD screen is easy to read and backlight included in replacement sensor models allow for taking measurements in almost any environment

Specifications

Accuracy:	2 mV
Backlight:	Yes
Battery requirements:	4 AAA
Casing IP rating:	IP67
Manual languages:	English, French, German, Italian, Spanish
Operating temperature range:	0 - 50 °C
Parameter:	ORP
Range:	±999 mV
Resolution:	1 mV
Sample volume:	9 mL
Warranty:	24 months
What's included?:	Includes: tester, 4 AAA batteries, lanyard and quick start guide.

What's included?

Includes: tester, 4 AAA batteries, lanyard and quick start guide.