



PHC3081-8 Combined pH electrode, Ag/AgCl, Temp.sensor, BNC

Product #: E16M305
AED Price: Contact Hach
Ships within 2 weeks

General-purpose pH electrode with Ag/AgCl reference system. Measuring range 0-12 pH, 0-80°C. Liquid junction: porous pin. With temperature sensor. BNC plug.

Specifications

Cable length:	1 m
Electrode junction :	Porous pin
Electrode type:	PHC3081-8
Filling solution:	Refillable: 3 M KCl + sat. Ag/Cl
Length:	103 mm
Material:	Glass
Measuring range:	0 - 12 pH
Sensor type :	Glass
Temperature range:	0 - 80 °C
Thermistor:	NTC 10 kohm
Warranty:	6 months

Required Accessories

- Filling Solution, Reference, 3 M KCl with AgCl, 100 mL (Item S21M004)