



Manual flow cell kit for 2100N/AN laboratory turbidimeters

| | |
|----------------------|----------------|
| Product #: | 4744900 |
| AED Price: | Contact Hach |
| Ships within 2 weeks | |

Low-pressure, flow-cell kit allows for a constant optical path.

Obtain colorimetric measurements with Hach spectrophotometers and turbidimeters without individual sample cells. Eliminate errors due to scratched or smudged sample cells. No need for optically matched sample cells. For rapid, convenient processing of large numbers of samples.

Specifications

| | |
|----------------|---|
| Includes: | Flow-cell stand assembly, flow-cell inlet reservoir, reservoir cover, collection drain assembly, flow-cell assembly, interconnecting tubing, flow-cell light shield |
| Sample volume: | 120 mL purge volume |
| Weight: | 1.9 kg |