



Automatic Cleaning Module for TU5300sc and TU5400sc

Product #: LQV159.98.00002
AED Price: Contact Hach

Cleaning your turbidimeter takes valuable time that your staff could use to complete other tasks. The Automatic Cleaning Module for the TU5300sc and TU5400sc turbidimeters completes cleaning for you, increasing instrument uptime.

Everything About Turbidity – Faster

The TU5 Series Automatic Cleaning Module can be set for hands-off cleaning of your sample cell.

Multiple Application Cleaning

Wipers are included with the cleaning module that can clean up to the most demanding finished water with iron and manganese.