



Maintenance Kit for TU5300sc and TU5400sc Laser Turbidimeter, with RFID Sealed Vials

Product #:

LZY907.98.00002 **Hazardous**

AED Price:

Contact Hach



Specifications

What's included?:

Calibration lid with user manual, lint-free cloths, wire maintenance rack, vial cleaning tool, vial compartment cleaning tool, replacement automatic cleaning module wiper, 20 NTU sealed vial Stablcal Standard with RFID, 600 NTU sealed vial Stablcal Standard with RFID, #0.1 NTU secondary Standard glass rod with RFID, vial removal tool

What's included?

Calibration lid with user manual, lint-free cloths, wire maintenance rack, vial cleaning tool, vial compartment cleaning tool, replacement automatic cleaning module wiper, 20 NTU sealed vial Stablcal Standard with RFID, 600 NTU sealed vial Stablcal Standard with RFID, #0.1 NTU secondary Standard glass rod with RFID, vial removal tool