



## Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

**Product #:**

**143601**

**Hazardous**

AED Price:  
Available

Contact Hach



Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as  $\text{CO}_2$ , with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

**Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.**

---

### Specifications

Case style:	D
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$ , ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Measuring range:	2 - 40 mg/L
Measuring range (2):	1.25 - 25 mg/L
Measuring range (3):	5 - 100 mg/L $\text{CO}_2$
Method:	Drop count titration/  Sodium hydroxide  Phenolphthalein
Number of tests:	200
Parameter:	Carbon Dioxide - As $\text{CO}_2$
Platform:	Drop Count Titration
Smallest increment:	1.25  2  5
Test kit model:	CA-23