



Chloride Low Range Test Kit, Model 8-P

Product #:

144001

AED Price:
Available

Contact Hach

Hazardous



Using silver nitrate titrant, this kit provides fast, accurate detection of chloride.

Low cost and easy to use drop count titration method using a silver nitrate titrant

Ideal for quickly determining quantitative results

Specifications

Case style:	B
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Measuring range:	5 - 100 mg/L & 20 - 400 mg/L Cl^-
Method:	Drop count titration/ Silver nitrate
Number of tests:	100
Parameter:	Chloride - As Cl^-
Platform:	Drop Count Titration
Smallest increment:	5
Test kit model:	8-P