



## TEST KIT, 5-EP mg/L HARDNESS

**Product #:**

**145401**

**Hazardous**

AED Price:  
Available

Contact Hach



Total Hardness drop count titration/EDTA test kit with range of 20 - 400 mg/L; contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Case style:	B
Measuring range:	20 - 400 mg/L $\text{CaCO}_3$
Method:	Drop count titration/  EDTA-powder pillows
Number of tests:	100
Parameter:	Hardness, total - As $\text{CaCO}_3$
Platform:	Drop Count Titration
Smallest increment:	20
Test kit model:	5-EP MG-L