



Total chlorine test kit (colour disk)

Product #:

1454202

Hazardous

AED Price:
Available

Contact Hach



Colour comparison test for easy and low-cost determination of total chlorine. Measuring range 0.02 - 0.7 / 0 - 3.4 mg/L Cl₂.

Specifications

Case style:	D
Measuring range:	Dual Range
Measuring range (2):	0.02 - 0.68 mg/L Cl ₂
Method:	Colour disc/DPD
Number of tests:	100 Total
Parameter:	Chlorine Total
Platform:	Color Disc
Smallest increment:	0.02
	0.1
Test kit model:	CN-70T