



## Bottle, BOD, 300 mL, 24/pk

<b>Product #:</b>	<b>62124</b>
AED Price:	Contact Hach
Available	

The standard for BOD analysis. Borosilicate clear glass bottles with flared lip and ground-glass stoppers to prevent air entrapment and for water seal. Frosted marking area. Stoppers included. Pack of 24.

---

### Specifications

Body Material:	Glass
Capacity:	300 mL
Material:	Borosilicate glass
Package size UOM:	24pk
Quantity:	24 /pk