



USB barcode hand-scanner for spectrophotometers

Product #:	LZV566
AED Price:	Contact Hach
Ships within 2 weeks	

The scanner is as an optional part of spectrophotometers for automatic identification (ProID) of norm and individual barcodes.

Recognition of different coloured cuvette-clips for the sample identification system ProID.

Specifications

Dimensions (H x W x D):	80 mm x 100 mm x 175 mm
Enclosure waterproof rating:	IP30
Interface:	1 USB
Max. operating humidity:	90 %
Operating temperature range:	0 - 55 °C
Quantity:	1 pcs
Storage conditions:	-20 °C - 70 °C
Warranty:	24 months
Weight:	0.185 kg