



Sension+ ORP/Redox liquid electrode, refillable

Product #: LZW5057.97.0002

AED Price: Contact Hach
Available

ORP/Redox electrode with liquid electrolyte and platinum wire. Glass body. Range ± 2000 mV, 0-80 °C. For Benchtop meter.

Specifications

Accuracy:	0.2mV with MM374 Meter
Cable length:	1 m
Connector:	BNC Connector
Dimensions (D x L):	12 mm x 130 mm
Electrode junction :	Ceramic pin
Electrode type:	Laboratory
Filling solution:	LZW9500.99
Length:	130 mm
Material:	Sensor Body: Glass
Measuring range:	± 2000 mV
Parameter:	ORP
Sensor type :	Platinum
Temperature range:	0 - 80 °C
Type:	Lab
Warranty:	12 months