



ORP gold electrode 5354, gel electrolyte, 0-80 °C, 6 bar, S8 screw head

Product #: LZU5354.97.0002
AED Price: Contact Hach

ORP gold electrode 5354 with gel electrolyte and S8 screw head. ± 2000 mV, 0 - 80 °C, up to 6 bar. Ceramic diaphragm.

Suitable for acid samples

containing Fe or Cr.

Specifications

Cable connection:	S8
Diameter:	12 mm
Diaphragm:	ceramic
Electrolyte:	gel
Length:	120 mm
Material:	Sensor: glass/Au
Measuring range:	± 2000 mV
Parameter:	ORP
Pressure range:	Up to 6 bar
Reference system:	Encapsulated Ag/AgCl crystals
Temperature range:	0 - 80 °C
Warranty:	6 months