



pH electrode 5338, plastic, gel electrolyte, fixed cable, 0-80 °C, 2 bar

Product #: LZU5338.97.0002
AED Price: Contact Hach

Low cost pH electrode 5338 with plastic body, gel electrolyte and fixed cable. 0 - 14 pH, 0 - 80 °C, up to 2 bar.

Specifications

Electrolyte:	Gel
Material:	Plastic
Measuring range pH:	0 - 14 pH
Parameter:	pH
Pressure range:	Up to 2 bar
Temperature range:	0 - 80 °C
Warranty:	6 months