



ORP platinum electrode 5358, plastic, gel electrolyte, fixed cable, 0-80 °C, 2 bar

Product #: LZU5358.97.0002
AED Price: Contact Hach

Low cost ORP platinum electrode 5358 with gel electrolyte and fixed cable. ± 2000 mV, 0 - 80 °C, up to 2 bar. Ceramic diaphragm.

Suitable for swimming pools,

drinking water, and aqueous samples in general.

Specifications

Cable connection:	Fix cable
Diaphragm:	ceramic
Electrolyte:	gel
Length:	120 mm
Material:	Sensor: external PC/internal glass/Pt
Measuring range:	± 2000 mV
Pressure range:	Up to 2 bar
Reference system:	Ag/AgCl wire
Temperature range:	0 - 80 °C
Warranty:	6 months