



ORP platinum electrode 5362, gel electrolyte, 3 diaphragms, 0-80 °C, 6 bar, S8 screw head

Product #: LZU5362.97.0002
AED Price: Contact Hach

ORP platinum electrode 5362 with gel electrolyte and S8 screw head. ± 2000 mV, 0 - 80 °C, up to 6 bar. 3 ceramic diaphragms.

Suitable for swimming pools,

drinking water, gas scrubbing towers, boiler water, detoxification, wwtp outlets etc.

Specifications

Cable connection:	S8
Diaphragm:	3 ceramic
Electrolyte:	gel
Length:	120 mm
Material:	Sensor: glass/Pt
Measuring range:	± 2000 mV
Pressure range:	Up to 6 bar
Reference system:	Encapsulated Ag/AgCl crystals
Temperature range:	0 - 80 °C
Warranty:	6 months