



## Conductivity cell 5388 with PT1000, k=0.6, 80 °C, MP-5, PG 13.5

**Product #:** LZU5388.97.0002  
**AED Price:** Contact Hach

Conductivity cell 5388 + PT1000 A.T.C., titanium, k=0.6, 10...80000 µS/cm, 0...80 °C. MP-5 connector, PG 13.5 thread.

---

### Specifications

Cable connection:	MP-5
Diameter:	12 mm
Length:	110 mm
Material:	Measuring Cell: Titania
Measuring range:	10 µS/cm - 80 mS/cm
Temperature compensation:	0 - 80 °C
Temperature range:	Pt 1000
Warranty:	6 months
What's included?:	Includes: sensor with cable and manual

---

### What's included?

Includes: sensor with cable and manual