



Conductivity cell 5390 with Pt1000, k=0.1, 100°C, MP-5, 1/2"

Product #: LZU5390.97.0002
AED Price: Contact Hach

Conductivity cell 5390 + PT1000 A.T.C., for pipe, in titanium, PVDF and stainless steel, 0,1...10000 $\mu\text{S}/\text{cm}$, 0...100 °C. MP-5 connector, 1/2" thread.

Specifications

Cable connection:	MP-5
Cell constant:	0.1 cm^{-1}
Length:	50 mm
Material:	Measuring Cell: Titanium
Measuring range:	0.10 - 10.000 mg/L ($\mu\text{S}/\text{cm}$)
Temperature compensation:	0 - 100 °C
Temperature range:	Pt 1000