



Conductivity cell 5395, k=0.1, 80 °C, S8 connector, PG 13.5

Product #: LZU5395.97.0002
AED Price: Contact Hach

Conductivity cell 5395, in glass and platinum, k=0.1, 0,1...50000 $\mu\text{S}/\text{cm}$, 0...80 °C. S8 connector, PG 13.5 thread. (Without A.T.C.).

Specifications

Cable connection:	S8
Cell constant:	0.1 cm^{-1}
Length:	110 mm
Material:	Measuring Cell: Glass/Pt
Measuring range:	0.10 - 50.000 mg/L ($\mu\text{S}/\text{cm}$)
Temperature compensation:	0 - 80 °C