



Conductivity cell 5396, k=1, 80 °C, S8 connector, PG 13.5

Product #: LZU5396.97.0002
AED Price: Contact Hach

Conductivity cell 5396, in glass and platinum, k=1, 1...200000 µS/cm, 0...80 °C. S8 connector, PG 13.5 thread. (Without A.T.C.)

Specifications

Cable connection:	S8
Cell constant:	1 cm ⁻¹
Length:	110 mm
Material:	Measuring Cell: Glass/Pt
Measuring range:	1 - 200.000 mg/L (µS/cm)
Temperature compensation:	0 - 80 °C