



Conductivity cell 5398 with PT1000, k=0.3, 80 °C, MP-5, PG 13.5

Product #: LZU5398.97.0002
AED Price: Contact Hach

Conductivity cell 5398 + PT1000 A.T.C., titanium, k=0.3, 0.5...80000 µS/cm, 0...80 °C. MP-5 connector, PG 13.5 thread.

Specifications

Cable connection:	MP-5
Diameter:	12 mm
Length:	110 mm
Material:	Measuring Cell: Titanio
Measuring range:	0.50 - 80.000 mg/L (µS/cm)
Temperature compensation:	0 - 80 °C
Temperature range:	Pt 1000
Warranty:	6 months
What's included?:	Includes: sensor with cable and manual

What's included?

Includes: sensor with cable and manual