



Conductivity cell 5399 with PT1000, k=10, 80 °C, MP-5, PG 13.5

Product #: LZU5399.97.0002
AED Price: Contact Hach

Conductivity cell 5399 + PT1000 A.T.C., glass and platinum, k=10, 500...500000 $\mu\text{S}/\text{cm}$, 0...80 °C. MP-5 connector, PG 13.5 thread.

Specifications

Cable connection:	MP-5
Cell constant:	10 cm^{-1}
Length:	110 mm
Material:	Measuring Cell: Glass/Pt
Measuring range:	100 - 500,000 mg/L ($\mu\text{S}/\text{cm}$)
Temperature compensation:	0 - 80 °C
Temperature range:	PT1000