



Ferrous Iron Reagent Powder Pillows, 0.02-3.00mg/L Fe (II)

Product #:

103769

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed packages for the determination of Iron in the ferrous form (Fe^{2+}). pk/100.

Specifications

Instrument:	All except PC II
Measuring range:	0.02 - 3.00 mg/L Fe^{2+}
Method:	1,10 Phenanthroline
Method number:	8146
Number of tests:	100
Parameter:	Iron, ferrous
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)