



Iron Reagent Powder Pillows (TPTZ), 0.012-1.800mg/L Fe

Product #:

2275699

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of low levels of iron. For 25mL sample. 100/pack

Specifications

| | |
|------------------|-----------------------|
| Instrument: | DR/3000, DR/2000 |
| Measuring range: | 0.012 - 1.800 mg/L Fe |
| Method: | TPTZ |
| Method number: | 8112 |
| Number of tests: | 100 |
| Parameter: | Iron |
| Platform: | Powder Pillows |
| Quantity: | 100 /pk |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)