



Iron Reagent, Accuvac, 0.012-1.800 mg/L Fe

Product #: 2510025

Hazardous

OBSOLETE ITEM
This item is no longer available.



Suggested replacements:

- 2608799

Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of low levels of iron by the TPTZ method.

Specifications

Instrument:	All except DR/820
Measuring range:	0.012 - 1.800 mg/L Fe
Method:	TPTZ
Method number:	8112
Number of tests:	25
Parameter:	Iron
Platform:	Accuvac
Quantity:	25 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)