



Molybdenum 2 Reagent Solution (100mL MDB)

Product #:

2352532

Hazardous

AED Price:
Available

Contact Hach



For the determination of Molybdenum, 0.02-3.00mg/L Mo. For 50mL sample. Measured dropping bottle.

Specifications

Description:	Replacement reagent for 2440200
Method:	Ternary Complex
Number of tests:	Varies
Platform:	Powder Pillows
Quantity:	100 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)