



Molybdenum Reagent Set, 20 mL, LR Molybdenum Reagent Set, 0.02-3.00 mg/L Mo

Product #:

2449400

Hazardous

AED Price:
Available

Contact Hach



For the determination of Molybdenum, ternary complex method. For 20mL sample.

Specifications

Instrument:	All except DR/3000, DR/2000, DR/820
Measuring range:	0.02 - 3.00 mg/L Mo
Method:	Ternary Complex
Method number:	8169
Number of tests:	100
Parameter:	Molybdenum, Molybdate
Sample size:	20 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)