



ChromaVer Chromium Reagent Powder Pillows, 0-0.60 mg/L

Cr⁶⁺

Product #:

1206699

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Chromium, hexavalent and total. For 25mL sample. 100/pack

Specifications

EPA compliant:	Yes
Instrument:	DR/3000, DR/2000
Measuring range:	0 - 1.00 mg/L Cr ⁶⁺
Method:	1,5-Diphenylcarbohydrazide
Method number:	8023
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)