



ZincoVer 5 Zinc Reagent Powder Pillows

Product #:

1403268

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for zinc determination. For 50 mL sample.

Specifications

EPA compliant:	Yes
Measuring range:	0.01 - 2.00 mg/L Zn
Method:	Zincon
Method number:	8009
Number of tests:	25
Parameter:	Zinc
Platform:	Powder Pillows
Quantity:	25/pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)