



CuVer 1 Copper Reagent Powder Pillows, 0.04-5.00mg/L Cu

Product #:

1418899

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Copper, Bicinchoninate Method. For 25mL sample. 100/pack

Specifications

EPA compliant:	Yes
Instrument:	DR/3000, DR/2000
Measuring range:	0.04 - 5.00 mg/L Cu
Method:	Bicinchoninate
Method number:	8506
Number of tests:	100
Parameter:	Copper
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)