



ECR Reagent Powder Pillows, pk/25

Product #:	2380268
AED Price:	Contact Hach
Available	

Powder reagent in sealed package for the determination of Aluminium, Eriochrome Cyanine R Method. For 25mL sample.

Specifications

Method:	Eriochrome Cyanine R
Method number:	8326
Method Number 1 Chemistry:	8326
Number of tests:	25
Parameter:	Aluminum
Platform:	Powder Pillows
Quantity:	25/pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)