



Manganese Reagent Set, Powder Pillows, periondate oxidation, 0.1-20.0 mg/L Mn, 10 mL sample

Product #:

2430000

Hazardous

AED Price:
Available

Contact Hach



Powder reagents in sealed packages for the determination of soluble Manganese, Periodate Oxidation Method. For 10mL sample. 100 tests

Specifications

EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000
Measuring range:	0.1 - 20.0 mg/L Mn
Method:	Periodate Oxidation
Method number:	8034
Number of tests:	100
Parameter:	Manganese
Platform:	Powder Pillows
Sample volume:	10 mL