



## Copper Reagent, Accuvac, 0.04-5.00 mg/L Cu

Product #: 2504025

Hazardous

**OBSOLETE ITEM**  
This item is no longer available.



Suggested replacements:

- 2188299

Ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Copper in presence of high levels of Aluminium, Calcium, and Iron.

---

### Specifications

Instrument:	All except DR/850, DR/820
Measuring range:	0.04 - 5.00 mg/L Cu
Method:	Bicinchoninate
Method number:	8026
Number of tests:	25
Parameter:	Copper
Platform:	Accuvac
Quantity:	25 /pk

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)