



Copper Masking Reagent Powder Pillows, 10 mL, pk/100

| | |
|-------------------|----------------|
| Product #: | 2603449 |
| AED Price: | Contact Hach |
| Available | |

Powder reagent in sealed package for the determination of Copper, Porphyrin Method. For 10ml sample.

For Copper determination by the Porphyrin method. Hach Method 8143 Range: 2 - 210 µg/L Cu. Sample Size: 10 mL Pack of 100 Powder Pillows.

Specifications

| | |
|------------------|----------------|
| Method: | Porphyrin |
| Number of tests: | 100 |
| Platform: | Powder Pillows |
| Quantity: | /pk |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)