



Potassium hydroxide standard solution, 8.00 N, 500 mL

Product #:

28249

Hazardous

AED Price:
Available

Contact Hach



For the determination of Copper, 0.1-8.0 mg/L Cu, and Hardness, Titration and Buret Titration Method.

Specifications

Concentration:	8.00 N
Method:	EDTA Titration
Parameter:	Varies
Platform:	Solution
Volume:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)